

# EuroSplint®

to immobilise broken bones  
and sprain injuries



Available in two sizes.  
Lightweight, washproof, compact,  
reusable and can be X-rayed.

## EuroSplint

The EuroSplint is an incredibly versatile product and is used to immobilise broken bones and sprain injuries. The splint can be shaped to suit the injured part of the body and held in place by loosely wrapping a bandage or tape around the splint. It can be used on the neck, shoulders, arms, fingers, wrist, knee, legs, and ankles. Special metal, macromolecule foam keeps the splint clean and dry. Store in well ventilated area at room temperature, away from direct sunlight. Can be used with X-Rays.

REF	PRODUCT	SIZE	UNIT	PACK
748	EuroSplint Roll	110mmL x 90mmW x 6mmD	1	1
2749	EuroSplint Folded	110mmL x 90mmW x 6mmD	1	1
2748	EuroSplint Finger Flat	95mmL x 45mmW x 6mmD	1	10

## EuroSplint - Finger

Used to immobilise injured fingers. The splint can be shaped around or on the finger, loosely wrapping a bandage or tape around the splint to hold in place. Washproof and reusable. Made from special metal, macromolecule foam.

REF	PRODUCT	SIZE	UNIT	PACK
2748	EuroSplint Finger Flat	95mmL x 45mmW x 6mmD	1	10

